

Poster Presentation

Open Access

Increasing the ease and speed of eukaryotic protein expression: a new cell-free in vitro translation system based on Sf insect cell extracts

Frank Schaefer*, Annette Zacharias, Nicole Brinker-Krieger and Uritza von Groll

Address: QIAGEN GmbH, QIAGEN Straße 1, 40724 Hilden, Germany

* Corresponding author

from The 4th Recombinant Protein Production Meeting: a comparative view on host physiology
Barcelona, Spain. 21–23 September 2006

Published: 10 October 2006

Microbial Cell Factories 2006, **5**(Suppl 1):P33 doi:10.1186/1475-2859-5-S1-P33

© 2006 Schaefer et al; licensee BioMed Central Ltd.

Background

For researchers looking for fast access to their protein of interest, cell-free expression systems are an attractive option, as they do not require specialized equipment, are open systems, and offer fast screening of expression constructs for expression efficiency, yield, solubility and other criteria (e.g. requirement of cofactors).

A broad range of eukaryotic proteins require posttranslational modifications such as phosphorylation, glycosylation, or signal peptide cleavage to display full functional activity.

Results

We will describe the most important features of a new system for recombinant eukaryotic protein expression that retains the speed and economy inherent to cell-free methods. Data showing the ability of the system to generate complex posttranslational modifications (PTMs), including highly efficient glycosylation of EPO and ORM1 pro-




Figure 1

6xHis-tagged EPO (erythropoietin) and ORM1 (Alpha-1-acid glycoprotein I) were expressed as ^{14}C -Leu-labeled proteins using the EasyXpress Protein Synthesis Insect Kit. After expression an aliquot of each was treated with the endoglycosidase EndoH, which removes glycan moieties from glycosylated proteins. Aliquots of treated and untreated protein separated by SDS-PAGE and proteins visualized by autoradiography.




Figure 2

Akt1 was expressed in triplicate reactions using the EasyXpress Protein Synthesis Insect Kit. 2.5 μl aliquots were separated by SDS-PAGE and protein detected by immunoblotting using an antibody specific for akt1 phosphorylated at Ser 473. NTC: no template control; M: marker.

teins (Figure 1) and phosphorylation of AKt1 kinase (Figure 2) will be presented. The robustness of the system is demonstrated by the expression of membrane proteins (e.g., OGCP) and the synthesis of several proteins, such as human clotting factors, that could not be expressed in *E. coli* systems.

Conclusion

Expressed proteins with the EasyXpress Protein Synthesis Insect Kit are post-translationally modified with high efficiency.

Publish with **BioMed Central** and every scientist can read your work free of charge

"BioMed Central will be the most significant development for disseminating the results of biomedical research in our lifetime."

Sir Paul Nurse, Cancer Research UK

Your research papers will be:

- available free of charge to the entire biomedical community
- peer reviewed and published immediately upon acceptance
- cited in PubMed and archived on PubMed Central
- yours — you keep the copyright

Submit your manuscript here:
http://www.biomedcentral.com/info/publishing_adv.asp

